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Rollator vital 54

Art.-Nr. 11448 001 grau
Art.-Nr. 11448 011 mokka

Rollator vital 62

Art.-Nr. 11448 000 grau
Art.-Nr. 11448 010 mokka

HMV-Nr. 10.50.04.1174

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Rollator vital 62

Art.-Nr. 11448 000 grey
Art.-Nr. 11448 010 mokka

Please read these notes carefully before use. Be sure to follow the instructions and keep them in a safe place! If this product is to be used by other persons, they must be fully informed of these safety warnings before use.

A rollator helps those with limited mobility maintain their existing mobility and enables them to enjoy a larger sphere of activities. Staying in the move with a rollator not only serves to improve circulation, coordination and balance, but it also enhances the user's well-being.

The Rollator was developed for use indoors and outdoors. It features variable and simple height adjustment of the handles, extremely easy brakes that function safely and puncture-resistant tyres with TPE coating. The frame and net seat are suitable for a maximum permissible loading of up to 150 kg. The Rollator can be folded for easy transport. It stands securely when folded.


Safety Indications

- Talk to your doctor or specialist retailer to see whether this Rollator is optimally suitable for you and your needs.
- Check that the Rollator is in perfect condition every time before you use it. Check the following components and make sure:
 - all screws, nuts and other parts are firmly tightened
 - all wheel surfaces are intact
 - the height adjustment screws are firmly tightened
 - the brakes function both as parking brakes and hand brakes
 - there are no visible signs of wear on the frames, welded joints and plastic elements.
 Do not use a damaged or defective Rollator!
- Ensure that the Rollator is correctly engaged before every use.
- Always activate the parking brake when you stop or park the Rollator. Always use the hand brakes slowly, evenly and both together.
- The braking effect can be reduced if the tyre surface is worn. Inspect the braking action of the brakes regularly.
- The Rollator should only be used on firm and even surfaces. There is a risk of slipping on dirty, wet, frosty or snow-covered ground! The braking effect is also limited. Please adapt the way you move and walk to these dangers.
- The Rollator should not be used as a support when getting up from a crouching or sitting position! It should only be used for its intended purpose and should not be used for transporting persons or objects over 5 kg.
- Ensure that the handles are adjusted to the correct height for the user in order to avoid improper weight load or falls. We recommend you walk between the handles. Always hold both handles; otherwise you can lose your balance.
- Make sure that your hands or clothes do not come into contact with the moving parts of the Rollator! Be particularly careful when folding and unfolding the Rollator.
- The Rollator should be used on the pavement. The Rollator should not be used in moving traffic. For your personal safety on the roads it is recommended that you also wear your own reflectors as well as those fixed to the Rollator.
- Be careful near children. Children should not play with the Rollator!
- The maximum permissible loading for the Rollator is 150 kg. The bag should only be loaded with up to 5 kg.
- Do not hang any handbags or carrier bags on the handles, the balance of the Rollator can be impaired if you do.
- Do not carry out any adjustments or modifications, which impair the safety of the design.
- Only use accessories approved by the manufacturer.
- Do not carry or pull the Rollator either by the strap to fold it up or by the brake cables.

Contents/Description of the parts:

- | | | | |
|----------------------------|----------------|-----------------------------------|------------------|
| A. Handle | E. Brake cable | I. Strap for folding the Rollator | K. Bag |
| B. Break lever | F. Side rails | J. Brake pads and edge guide | L. Bag holder |
| C. Seat | G. Frame | | M. Magnetic lock |
| D. Height adjustment screw | H. Step aid | | |


Assembly

Remove all packaging. The packaging should be disposed of responsibly. The plastic and cardboard parts of the packaging can be recycled.


Opening the Rollator

1. Pull the handles apart.
2. Push the side rails downwards until they engage. The Rollator is now fully opened.
3. Make sure the Rollator is engaged before you use it.

Folding up the Rollator


1. Fold the Rollator up by pulling the strap upwards.
2. Push the handles inwards and the Rollator together. The fully folded Rollator is locked with magnets.

Height adjustment of the handles

Adjust the height of the handles so that the user can walk upright, without bending over. It is recommended that the user stands upright with arms at their sides between the handles. The grips should be adjusted to the height of their wrists.


1. Loosen the height adjustment screw A.
2. Move the handle to the desired height.
3. Tighten the height adjustment screw A.

Brake

Check the brakes function before use.


1. Speed is reduced by pulling the brake lever upwards.
2. The parking brake is activated by pressing the brake lever down until it engages.
3. The parking brake is released from its locked position by pulling the brake lever upwards.
4. The Rollator should not be pushed along with the parking brake on. If the Rollator is pushed along with the parking brake on this can cause extreme wear to the tyres and brakes and their function impaired.
5. The Rollator must not be used with damaged brakes.

Brake adjustment

The brakes must be adjusted so that both wheels brake evenly. The distance between the brake pad and the wheel surface should be about 2 mm.


1. Loosen the locking nut.
2. Adjust the braking force by turning the brake cable spanner inwards or outwards.
3. Check the braking function by turning the wheels and operating the brake levers. The wheels should turn freely but stop immediately when the brakes are applied.
4. Tighten the locking nut firmly again.

Seat


Activate the parking brake before you sit on the seat. Always sit fully on the seat with your back in the direction of travel. Do not lean against the seat! That could impair the stability of the Rollator.

The seat has been tested for a user weight of up to 150 kg.

Negotiating obstacles with the support of the step aid

Tread on the step aid in order to raise the front wheels of the Rollator so that you can negotiate a doorstep, kerb or another obstacle more easily.


Care and Maintenance

Rollator

The Rollator may only be used as a walking aid. Keep the Rollator free from dirt, so that it remains functional. The Rollator can be cleaned with a neutral cleaning agent and a cloth or brush and then rinsed and dried. A dirty Rollator can also be cleaned with a high pressure cleaner or hosepipe. However, you should avoid spraying water directly into the ball bearings. Dry the Rollator thoroughly after cleaning or exposure to rain. The brake pads and springs should also be cleaned regularly.

The Rollator should only be disinfected by a specialist. The product can be disinfected with a disinfectant with 70-80 percent ethanol content. Disinfectants containing chlorine or phenol should not be used.

Do not use the Rollator in temperatures below -35°C or above $+60^{\circ}\text{C}$. Apart from small corrections on the paintwork, changing wheels or adjusting the brakes all repairs should be carried out by a specialist retailer.

Bag

Re-use

The product can be re-used. It must be cleaned and disinfected according to instructions. It must be inspected by specialist shops for any faults and then authorised before it can be used again. Please order our hygiene and maintenance guidelines.

Storage

Store the Rollator light and dust-protected in a dry place. ☂️🌞

Disposal

Dispose of this product with your local disposal enterprise.

Technical Details

Material:

- Frame: Aluminium with plastic-segments
- Wheels: TPE

Rollator vital 62

Dimensions:

- 60 x 66 x 74-101,5 cm (W x L x H)
- Folding measure: 23 x 66 x 79
- Wheels: 20,2 x 3,5 cm
- Seat level: 62 cm
- Distance between handles: 44 cm
- Weight: 6,8 kg + bag (300 g)
- Turning circle: 84,5 cm
- Recommended body size: 130-210 cm

Bag

- Material: 100 % Polyester
- Dimensions: 32 x 15 x 28 cm (W x L x H)
- Permissible weight: max. 5 kg

Rollator vital 54

Dimensions:

- 60 x 66 x 66-85 cm (W x L x H)
- Folding measure: 23 x 66 x 73,5 cm
- Wheels: 20,2 x 3,5 cm
- Seat level: 54 cm
- Distance between handles: 44 cm
- Weight: 6,64 kg + bag (300 g)

CE Mark

The Rollator bears the CE mark and complies with EU Council Directive 93/42/EEC concerning medical devices. The type label can be found on the frame and contains the following information:

- CE mark
- Permissible weight
- Serial number
- Manufacturer details

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