Bedienungsanleitung | Instruction manual  
Manual de instrucciones | Mode d'emploi

## Centuro S7


DE

EN

ES

FR

*Enjoy mobility.*

## Inhalt

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## Table of contents

1.	Preface	38	4.	Using the scooter	50
1.1.	Signs and symbols	38	4.1.	Control Functions	50
2.	Important safety instructions	39	4.1.1.	Changes to the controller	53
2.1.	General safety instructions	39	4.2.	Getting in and out	53
2.2.	Information on EMC interference	41	4.3.	Driving the scooter	54
2.3.	Brakes	41	4.4.	Pushing mode	55
2.4.	Safety in vehicles	41	4.4.1.	Brakes in pushing mode	56
2.5.	Use on Public Roads	42	4.4.2.	Initiating emergency stop	56
2.6.	Purpose	42	4.5.	Overload switch	56
2.7.	Indications	42	4.6.	Preparing the scooter for transportation	58
2.8.	Contraindications	42	4.7.	Disassembling for storage	59
2.9.	Declaration of conformity	43	4.8.	Taking on public transport	59
2.10.	Information about range	43	5.	Technical Data	60
3.	Description of Product and Contents	44	6.	Care instructions for the user	62
3.1.	Inspection on Delivery	44	6.1.	Cleaning and disinfection	62
3.2.	Type plate and serial number	45	6.2.	Check Before Use	63
3.3.	Contents	46	6.2.1.	Brake Inspection	63
3.4.	Summary	46	6.2.2.	Checking the lights, indicators and horn	63
3.5.	Basic adjustments	47	6.2.3.	Inspection of Tyres	63
3.5.1.	Steering column	47	6.3.	Batteries	64
3.5.2.	Seat	47	6.3.1.	Charging	64
3.5.3.	XL seat (optional)	47	6.3.2.	Replacing the batteries	65
3.5.4.	Back rest	48	7.	Maintenance tips for retailers	66
3.5.5.	Arm rest	49	8.	Storage / Shipping	66
3.5.6.	Seat belt (optional)	49	9.	Disposal	67
			10.	Use by others	67
			11.	Maintenance	67
			12.	Service life	68
			13.	Warranty	68

## 1. Preface

Dear user,

You have decided on a high-quality Bischoff & Bischoff scooter, and we would like to thank you for your trust.

The scooter was designed to offer you many advantages and meet your needs. Thanks to the various adjustment options and optimum adaptability, the scooter can be adjusted according to your individual body size. The scooter is designed for outdoor use only.

Please read and follow the operating manual carefully before using your new scooter for the first time. The operating manual is an integral and necessary component of the scooter. Keep the operating manual where you can refer to it and if you transfer the scooter, give the manual to the new owner.

The Scooter S7 is available with a max. speed of 15 km/h.

The scooter complies with the requirements of ISO 7176-14: 2008 (drive and controller systems for scooters and scooters; requirements and testing). For visually impaired users, this document is available as a PDF file at our website [www.bischoff-bischoff.com](http://www.bischoff-bischoff.com).

Any repairs and adjustments require special technical training and may therefore only be carried out by authorized Bischoff & Bischoff dealers.

### 1.1. Signs and symbols

Important! This is especially useful information on a given subject.

Warning! These are especially safety-relevant instructions. Follow the directions in the manual!

## 2. Important safety instructions

### 2.1. General safety instructions

- To avoid falls and dangerous situations, you should first practice driving your new scooter on level, unobstructed ground. We recommend having an assistant with you.
- Before operating the scooter, please check that all attached parts are correctly fastened.
- Practice driving on level, unobstructed ground before first use. Get a good understanding for the way the scooter brakes and accelerates when driving straight ahead or a curve.
- Please note that displacing the center of gravity by shifting the body or by adding loads may increase the risk of tipping over the scooter to the side or backwards.
- The scooter may not be used for transport of more than one person or any additional load. Make sure to only use it for the intended purpose.
- Be aware of the risk of trapping body parts when making settings and adjustments to the scooter.
- Never use the scooter under the influence of alcohol or other substances that may impair your attention and/or your physical and mental capacity.
- To overcome obstacles, please be sure to use ramps. Should it be necessary to lift the scooter over an obstacle, then only do so without the occupant. Only lift by the solid parts of the frame.
- Avoid hitting obstacles (steps, curbs) at full speed.
- The control unit must always be switched off when getting in or out of the scooter.
- All road traffic regulations must be observed when driving on public roads.

- The rear lights at the scooter must always be visible to other road users.
- If possible, wear bright, conspicuous clothing. This will make it easier for other road users to see you.
- Please note that the upholstery will heat up if directly exposed to sunlight and may cause injuries to your skin when touched. You should therefore cover these parts or protect the scooter from exposure to direct sunlight.
- The seat materials comply with the flammability resistance requirements of EN 1021-1 and EN 1021-2. Nevertheless, you should exercise the necessary caution if you are smoking and refrain from putting any cigarettes down on the scooter.
- Reporting incidents. The operator or user must immediately report
  1. any malfunction
  2. any change to the features or performance as well as
  3. Immediately report any improper labelling or instructions for use of a medical device that have or could have led to the death or serious deterioration of the health of a patient, an employee or a third party to your competent regional authority for medical devices. This institute will immediately forward the notification to the operator's competent authority and also inform the manufacturer and the manufacturer's competent authority.

Please note that operation in pushing mode by disconnecting the drives must only take place on a level surface. The electro-mechanical brake is thereby disabled. There is a danger of fatal injury if this instruction is not complied with..


## 2.2. Information on EMC interference

Despite compliance with all applicable EMC directives and standards, it is possible that the scooter may affect or be affected by other electrical devices (e.g. mobile phones, alarm systems in department stores).

If you become aware of any such interference, either turn off the other device or turn off your scooter and move it outside the range of the interference.

## 2.3. Brakes

There is a service brake for braking while in pushing mode with the drive disconnected or in an emergency (Fig. 01). After letting go of the joy-stick, the automatic brakes will slow the scooter until it comes to a complete stop. If the drive is connected, the scooter cannot be pushed or displaced because the brake is automatically applied to the motors.


The braking distance is significantly longer on a downhill gradient than on level ground. Keep this in mind when driving on sloping terrain and brake in time.

## 2.4. Safety in vehicles

The Centuro S7 is not suitable for use as a car seat. Under no circumstances may it be used as such. It does not meet the requirements stipulated by ISO 7176-19. Failure to observe this instruction can lead to extremely serious or fatal injuries in the event of an accident.

Your scooter is not suitable for use on the designated buses in the public transport network.

## 2.5. Use on Public Roads

The scooter is designed for outdoor use. Please keep in mind that you are using public roads and therefore have to observe all road traffic regulations.

Do not endanger other road users through inconsiderate driving, especially on walkways.

## 2.6. Purpose

The scooter Centuro S7 is designed for outdoor use by persons with impaired mobility. It is therefore capable of travelling long distances and overcoming obstacles outdoors. The scooter Centuro S7 is classified according to EN 12184 in category C. The maximum load capacity is 160 kg.

## 2.7. Indications

Inability or severely impaired ability to walk due to

- paralysis
- loss of limb
- defect/deformity of limb
- joint contracture
- joint damage (not in both arms)
- other diseases

Provision of a scooter is advisable if the use of manually operated scooters is no longer possible because of the disability but the person can still competently operate a vehicle driven by an electric motor.

## 2.8. Contraindications

The provision of a scooter is unsuitable for persons

- with severe balance disorders
- with impaired or insufficient vision
- with strongly impaired cognitive abilities

## 2.9. Declaration of conformity

Bischoff & Bischoff GmbH as the manufacturer declares under its sole responsibility that the Centuro S7 scooter complies with the requirements of Directive 93/42/EEC.

## 2.10. Information about range

The range of the scooter Centuro S7 is approx. 45 km. The range that can be achieved by electric scooters depends on the following factors:

- Composition of the carriageway
- Way of driving
- Battery condition
- Vehicle load weight
- Ambient temperature
- The use of accessories

Please note that the range specified by us was determined under test conditions:

- Batteries as new and fully charged
- Ambient temperature of 21 °C
- max. user weight of 160 kg
- Level, firm surface
- Constant speed
- Deactivated lighting system

We wish to advise you that the following factors severely restrict the range:

- Poor state of charge of batteries
- Low ambient temperature (below 10 °C)
- Frequent starting and braking
- Driving uphill
- Driving on surfaces that are not firm (field or forest paths)
- Tyre pressure incorrect or too low (see chapter 5)

If there is a combination of several factors, the range can be reduced by up to 50 %.

Never use your scooter with flat or almost dead batteries. You could put yourself and others at great risk through an unforeseen halt.

### 3. Description of Product and Contents

#### 3.1. Inspection on Delivery

All B+B products undergo a proper final inspection in our company and are labeled with the CE mark.

The scooter is delivered in a special carton. Keep the carton after unpacking, if possible. It is useful for later storage of the product or for sending it back.

Delivery and first instructions are taken care of by your qualified medical supply dealer.


If the scooter is shipped by rail or road, the goods should be checked immediately for transport damages (carton) in the presence of the bearer.

Should there be any damage, please contact our customer service (Tel.: +49 700 6000 7070) immediately.

Check the contents for completeness and make sure nothing is damaged. Should you find any irregularities or damage, please contact our customer service (Tel.: +49 700 6000 7070).

### 3.2. Type plate and serial number

The type plate (Fig. 02) and serial number are located on the rod of the seat.

- A Manufacturer logo
- B Manufacturer information
- C Medical device
- D CE mark
- E RoHS-Logo
- F Caution! - Follow instructions for us
- G Important - Follow instructions for use
- H Not suitable for use as a seat in a vehicle
- I Model designation
- J Model number
- K Application class
- L Date of manufacture
- M Maximum user weight
- N Maximum speed
- O Maximum weight capacity
- P UDI number
- Q Serial number

### 3.3. Contents

After receipt of the goods, please check the contents for completeness right away. The contents consist of:

- outer package
- Scooter (ready for use incl. 2 batteries)
- Charger
- Operating manual
- Seat
- Armrest
- Headrest
- Basket

### 3.4. Summary

### 3.5. Basic adjustments

Before using the scooter for the first time, it is necessary to adjust the scooter according to your height. This is essential to ensure that your posture is physiologically correct for driving for extended periods of time without excessive fatigue. The basic adjustments of the scooter are described below. You should carry out these adjustments with the help of your medical supply dealer.

#### 3.5.1. Steering column

Adjusting the steering column of your scooter enables you to set the driving position so that it is comfortable. You can make the adjustment while seated in the driver's seat.

- Hold the handlebar with one hand and activate the adjusting lever with the other hand (Fig. 03).
- Move the steering column into the required position and release the lever.

#### 3.5.2. Seat

The seat height can be adjusted and the seat can also be moved forwards and backwards as well as rotated to the side.

##### **Rotating the seat**

- Push down the lever (located on the right underneath the seat surface, Fig. 4). The seat can now be rotated in both directions.
- There is a fixed lock-in position for the seat every 45°.
- Release the lever and make sure that the seat is locked in position again.

#### 3.5.3. XL seat (optional)

As an option, an XL seat is available for the Centuro S7 (see Chapter 5, Technical Data). A compatible seatbelt is also optionally available for the XL seat.

**Moving the seat**

- Pull the lever (on the right below the seat surface, Fig. 05) towards you. The seat can now be moved forwards and backwards.
- You can now adjust the seat to the required position.
- Release the lever and ensure that the seat is locked in one of the closest fixed positions.

**Seat height**

You will firstly need to remove the seat.

- Fold the back rest onto the seat by holding the back rest with your right hand and using your left hand to lift up the lever on the left next to the seat surface (Fig. 08). The back rest will then fold forwards automatically.
- Simply remove the seat by pulling it upwards (Fig. 6).
- Undo the screw connection (Fig. 7) and position the seatpost in the required position inside the tube guide.
- Then securely tighten the screw connection again (Fig. 7) using 40Nm in the next hole.
- Now place the seat back onto the mounting.
- Make sure that it is locked in position again.

**3.5.4.Back rest**

The angle of the back rest on your scooter can be adjusted forwards and backwards. You also have the possibility to limit the adjustment range backwards

**Adjusting the angle**

- Pull the lever (on the left on the seat, Fig. 08) upwards. The back now pushes forwards.
- Push the back rest into the required position.
- Release the lever and make sure that the back rest is locked in position again.
- Please note that the adjustment of the seat to the rear and/or the angle adjustment to the rear can increase the overall length of the scooter.


**3.5.5.Arm rest**

The arm rests can be swivelled upwards in order to make it easier to mount and dismount. In addition, the height, distance and angle of the arm rests can be adjusted.

**Height**

- Loosen the screw until you can remove it, hold the screw in its position. Now you can move the arm rest vertically up and down (Fig. 9).
- Move the arm rest to its preferred height and tighten the screw.

**Side distance**

- Undo the hand screw (Fig. 10) on the back of the back rest frame while holding onto the armrest.
- Adjust to the required distance and tighten the screw again.
- The armrests can also be removed (e.g. for replacement, storage or transportation) after undoing the screw on the back of the frame (Fig. 10).
- Please note that pulling the armrests out to the side can increase the vehicle width.

**Angle of support**

- You can adjust the angle of support of the armrests by screwing the knurled nut (Fig. 11) on the relevant armrest further in or by unscrewing it.

**3.5.6.Seat belt (optional)**

The seat belt will be attached to the frame from the side.

- Adjust the belt length to meet your needs by positioning the clasp accordingly.
- Fasten the belt with the belt buckle.
- Unfasten the belt by pressing the red button (Fig. 12) on the belt buckle



## 4. Using the scooter

After you have adjusted the scooter according to your height, familiarise yourself with the functions of the controller.

### 4.1. Control Functions

The controller is a component of the steering column of your scooter. It contains LED displays and functional elements.

#### Hazard warning lights

Use this button (Fig. 14) to switch on the hazard warning lights in emergency situations where it is not possible to drive the scooter. All four indicator lights then flash and signal to other road users that your scooter cannot be driven.

#### Speed setting


You can use these buttons (Fig. 15) to set the maximum speed.

- Turn the switch to the left (turtle) to reduce the maximum allowed speed.
- Turn the switch to the right (rabbit), to raise up to the maximum allowed speed.
- The H/L switch allows a default setting for slow (low-L) and/or fast (high-H). Set at low, the speed is approx. 50% of the max. speed and at high, a max. of 15 km/h.
- To switch between km/h and mph, press the MODE button until the speed indicator flashes. Use the H/L key to switch between units. To confirm your selection, hold down the SET button for at least 3 seconds.

#### Headlight


Use this button (Fig. 16) to switch on the lights on your scooter when dusk begins to fall or when passing through dark routes such as tunnels. This will make you more easily visible to other road users.

#### Battery Status Display

The battery status display (Fig. 17) shows to what extent the batteries are charged. If all ten bars are visible, the battery is fully charged.


As you use the scooter, the battery capacity decreases and fewer bars are displayed. If only two bars are displayed, this means that only short journeys can now be made and the battery should be recharged (see charging).

#### Indicators (left/right)

Use this button (Fig. 18) to switch on the left or right indicator on your scooter to signal to other road users that you wish to turn left/right. Pressing the button again turns the indicator off again. After flashing 15 times the signal turns off automatically.

**Horn**

The horn (Fig. 19) enables you to draw attention to yourself in difficult situations.

**Ignition key**


Switch the scooter off and on using the ignition key. (Fig. 20).

**USB Port**

There is a USB port on the inside of the steering column, below the handlebar (Fig. 21). You can use this to conveniently charge your smartphone or other mobile terminals when your scooter is switched on. The max. current supply is 1 A/5 V.

**Mileage**

You can display your mileage in kilometers (km) or in miles (mi). Press the MODE button twice until the mileage indicator flashes. Use the H/L key to switch between units. To confirm your selection, hold down the SET button for at least 3 seconds.

**Time**

To set the time, press the MODE button three times until the time flashes. You can switch between „AM/PM“ and a 24-hour display by pressing the H/L key.

Set the hours by pressing the left arrow key (blinker). Set the minutes using the right arrow key (blinker).

To confirm your selection, hold down the SET button for at least 3 seconds.

**Temperature**

You can choose to display the ambient temperature in degrees Celsius (°C) or Fahrenheit (F). Press the MODE button four times until the temperature indicator flashes. Use the H/L key to switch between °C and F. To confirm your selection, hold down the SET button for at least 3 seconds.

After pressing the MODE button, the selected display flashes three times. If you do not make any changes during this time, you will leave the menu automatically.

**4.1.1. Changes to the controller**

Changes to the programming of the controller may only be made by an authorized supplier. Changes to the programming have an effect on the speed, acceleration and deceleration. This can lead to life-threatening injuries of the user and / or third persons if improperly carried out

**4.2. Getting in and out**

Ensure that the scooter is switched off before sitting on it.

- Move the steering column, as described previously, into the upright position and swivel back the arm rest on the side from which you intend to sit down so that the arm rest is out of the way.
- It is easier for you to sit down if you rotate the seat 90° to the side from which you intend to sit down (see the section entitled Basic adjustments).
- Now sit down.
- Rotate the seat in the driving direction or swivel the arm rest forwards again.
- Adjust the steering column so that you can comfortably reach all the controls.
- To get out, follow the instructions accordingly in reverse order. However, take out the key beforehand so that you can ensure that no unintended actions or movements occur.

Ensure that the seat firmly locks in position again

### 4.3. Driving the scooter

You should make your first journey on a level site with ample space

- Ensure that you are sitting correctly on your scooter. Set the speed as low as possible. You can increase the speed when you are more familiar with handling the scooter.
- To switch on your scooter, insert the ignition key into the ignition lock on the control panel and turn the key a quarter turn to the right.
- To switch off the scooter, turn the key a quarter turn to the left. After getting out of the scooter, always ensure that you have removed the ignition key and do not leave the scooter unattended with the ignition key inserted.
- Please make sure that you do not have your hand on the forwards / backwards lever when you switch on your scooter. This would cause your scooter to display an error message. In this case, turn it off and then on again to rectify this error.
- Pull back the speed selection lever on the right to move forwards (Fig 22). The further you pull back the lever, the higher the speed. The speed can go up to the maximum speed that has been preset using the speed control. Release the lever to stop gently. Practice these two procedures in order to get used to them.
- Steering is easy and works exactly like a bicycle in that you turn the handlebars in the direction in which you wish to move. Leave plenty of distance when you travel around a curve so that the rear wheels are able to follow the path of the front wheels and so that you do not cut across objects.


It is necessary to pay attention when reversing (fig 23). Ensure that your speed is low before reversing. Please keep in mind that you need to steer in the opposite direction when you are reversing. This requires practice. Practice this also on a level site with ample space.

Please note that, for safety reasons, the reversing speed is only half the speed going forward.

- If you turn sharply, turn the handlebars before you accelerate. Always drive in a wide arc around corners and obstacles and ensure that you drive slowly and carefully. For your own safety, the scooter automatically reduces the speed when cornering.
- When the drive is connected, the brake is automatically applied when the throttle lever is released. The scooter then cannot be pushed or displaced because the brake is automatically applied to the motor.

### 4.4. Pushing mode

If you want to push the scooter in certain situations, you must unlock the brakes. It is easy to push the scooter in free wheel mode.

The electro-mechanical brake will then not operate.

- Pull the unlocking lever in the back (Fig 24) downwards.
- A small jerk will indicate that the drives are now unlocked.
- To lock, simply push the lever up and locking then takes place automatically.
- The engine has a centrifugal clutch, the brake blocks automatically should you accelerate the scooter too hard in overrun cutoff. Only after a coming to a complete standstill will it be released.

For safety reasons, the motor drive is locked when the scooter is being operated in pushing mode.

Caution! Danger of fatal injury! The drives must only be disconnected on a level surface.

#### 4.4.1.Brakes in pushing mode

A manual service brake is attached to your scooter (as described in section 2.3). This can be used for support with braking when the drives are connected; it can also be used as a service brake when the drive is disconnected. The service brake acts on both front wheels and also serves as a locking brake.

Parking and locking the brake is performed by locking the drive. For this purpose, the scooter must be brought to a standstill by a helper and/or by using the manual service brake. The electro-mechanical brake (fig. 26), must then be activated.

#### 4.4.2.Initiating emergency stop

Should the automatic brake not work as usual while driving, you can initiate an emergency stop. Proceed as follows:

Turn the key a quarter turn to the left to stop the engine. The rear wheels will block immediately. Note that you can be thrown forward by the abrupt stop. Hold on to the handlebars adequately and perform this emergency stop only in an absolute emergency.


#### 4.5.Overload switch

The controller will cut out if the motor is overloaded, e.g. by driving uphill for too long. In this case, the ignition just needs to be switched off and on again in order to start moving again. If this does not work,

this means that the controller has overheated. In this case, please wait a few minutes. In the event of an electrical overload, the fuse will shut off the power to the scooter. The overload switch (Fig. 27) is located on the front of the seat construction (Fig. 28). After pressing the switch, the scooter should function normally again. If it does not function as expected, please contact your specialist retailer or Bischoff & Bischoff customer service.

**RISK OF ENTRAPMENT!** Ensure that you are aware of the increased risk of entrapment when stowing your scooter.


#### 4.6. Preparing the scooter for transportation

You can fold up the scooter in order to save space during transportation.

- Switch off the scooter.
- -Remove the armrests (Fig. 29) and then fold the back rest forwards onto the seat as described in Section 3.5.2.
- Remove the basket (Fig. 30).
- Fold the steering column towards the foot area (Fig. 31).
- Depending on the available space, it is sufficient to fold down the back rest and steering column and secure the scooter using securing straps on the frame (Fig. 32).
- Carefully protect the scooter against slipping by attaching the securing straps at the designated holding points.
- The scooter, including the batteries, is generally suitable for transportation on an aircraft. You can find out the transportation conditions from your airline. The batteries may need to be removed. Information on this can be found in Section 6.3.

#### 4.7. Disassembling for storage

For storage, you must also remove the batteries, be extremely careful. Avoid handling metal tools, touching the battery terminals and observe the instructions of the battery manufacturer.

This work should only be done by a trained technician. There is danger of fatal injury!

For your safety, the seat must be lifted by two people, as it is the heaviest individual part.

- You will firstly need to remove the seat and the panel.
- Remove the seat by pulling it upwards.
- The panel is fixed in place with 2 hand screws.
- Undo these screws and remove the panel.
- Disconnect the plug connection for the rear lighting.
- The batteries should then be removed. To make this easier, the cables are fitted with plug connections.
- Disconnect the plug connections (Fig. 34) for both batteries.
- Then undo the screw on the mounting bracket (Fig. 35). You do not need to completely remove the screw in order to do this.

#### 4.8. Taking on public transport

Your scooter is not suitable to be taken on public transport and is not authorised to bear the corresponding sticker (Fig. A).

## 5. Technical Data

### Dimensions and weights

Max. load:	160kg
Seat width:	457mm (XL seat: 520mm)
Max. effective seat width:	Seat width + 100mm
Seat depth:	457mm (XL seat: 550mm)
Seat angle:	0°
Seat height, from ground:	665/690/715mm of arm rest to seat surface
Distance armrest to seat:	230/280mm
Back rest height:	482mm, 665mm with headrest (XL seat: 550mm, 730mm with headrest)
Back rest angle:	30°
Overall length:	1440mm
Overall width:	650mm
Packing size Length/Height:	1365/640mm
Unloaden weight:	150kg
Weight of heaviest parts:	124kg, incl. batteries
Seat:	16,6kg
Head rest:	0,6kg
Arm rests (2 pcs.):	4,4kg

### Electric

Fuses:	5A (reset-able)
Batteries:	2 x 12V, 80A (maintenance-free, leak-proof)
Charger:	8A
Motor:	1100 W
Peak power:	1950 W
Controller:	PG 140A (S-DRIVE)
USB ports:	ca. 1A/5V

### Tyres

Driving and steering wheels:	13" (D330 mm), 250kPa (2,5bar)
------------------------------	-----------------------------------

### Driving characteristics

Speed (forwards):	15 km/h
Climbing capacity:	12°
Max. obstacle height:	100mm
Range:	ca. 45km
Reversing width:	2300mm
Braking distance (from max. speed):	< 1600mm
Static stability:	25° uphill 9° downhill 15° sideways

### Operating forces

Throttle lever:	< 13,5N
Buttons on the controller:	< 5N
Disconnecting the drives:	< 60N

### Operating conditions

Temperature:	-25 - +50°C
Air humidity:	20 - 80%

### Storage conditions

Temperature:	-40 - +65°C
Air humidity:	45 - 60%

### Equipment

Frame:	Powder-coated steel frame
Panel:	Plastic
Seat:	Artificial leather

2 independent braking systems

## 6. Care instructions for the user

### 6.1. Cleaning and disinfection

#### **Cleaning**

- Make sure that the scooter is switched off before starting to clean it.
- The frame components and panels of the scooter can be wiped with a damp cloth. You can use a standard mild cleaning agent if they are very dirty.
- The wheels can be cleaned with a damp brush with plastic bristles (do not use a wire brush).
- You can use a damp cloth to remove dust and minor dirt from the seat, back rest and armrests. Do not use aggressive cleaning agents because they may attack the surface of the material. The upholstery can be washed with warm water at 40 °C.
- Make sure that the electrical and electronic components do not come into contact with water unnecessarily.

#### **Disinfection**

- The scooter must be thoroughly cleaned in accordance with the instructions prior to disinfection.
- Use a standard disinfectant for wipe disinfection.
- Spray the upholstery and handles.
- Allow the disinfectant sufficient time to work before wiping off, as specified in the manufacturer's instructions.
- Make sure that the electrical and electronic components do not come into contact with liquid unnecessarily.

**Warning!**  
Do not use a high-pressure or steam cleaner.

### 6.2. Check Before Use

#### 6.2.1. Brake Inspection

Make sure your brakes are operating properly before each use.

- Start off carefully and let go of the throttle lever. The scooter must come to a complete stop as usual.
- When stopped, the mechanical brakes must engage with an audible click. Should you notice any unusual braking behavior, please notify your dealer right away and ensure the scooter is not being used.
- The motor brake is functioning correctly if you can't push the scooter when it is turned off. Should you nevertheless be able to do so, then the electro-magnetic brakes are defective. In this case you may no longer use the scooter. Contact your dealer immediately.

#### 6.2.2. Checking the lights, indicators and horn

Make sure your lights, indicators and horn are functioning properly before each use.

- If a unit under inspection does not react when the corresponding button on the control panel is pushed, please notify your dealer immediately and ensure the scooter is not being used.

#### 6.2.3. Inspection of Tyres

Before embarking on any journey, check that the tyres are in good condition.

- Check the tyre pressure. This should be 2.5 bar.
- Carry out a regular check of the tread depth of the tyres once a month. This should be at least 1 mm.
- Any faulty or worn out tyres must be changed through your mobility dealer.

### 6.3. Batteries

Your scooter is equipped with maintenance-free gel batteries whose fluid level doesn't need to be checked.

If your scooter is not used for an extended period, the charger should be connected about every 4 weeks to prevent damage to the batteries by deep discharge.

Before embarking on any journey, just check the charge status of your battery on the battery status display located on your controller. If the battery status display is in the low range, do not travel longer distances and promptly connect the battery charger to the scooter. We recommend connecting the charger to the scooter after every journey

#### 6.3.1. Charging

- Only use approved chargers (of the type specified in the technical data or supplied by the B+B specialist retailer) to charge the batteries.

Make sure to also follow the charger's instructions for use!

- Before switching on the charger, check that the power supply is 230 V.
- First of all, switch the scooter off.
- Connect the charger to the power socket (230 V)
- Check the LED status on the charger. If a red LED illuminates or flashes, do not start the charging process.
- If the charger is ready for use, connect it to the scooter. The charging socket is located on the side underneath the ignition lock (Fig. 40).
- Upper LED is yellow: Charger is switched on  
Lower LED is flashing green: Charger ready for use.
  - Orange: Charging in process
  - Orange/green: Charging in progress, over 80% charged
  - Green: Charging complete
  - Flashing red: Error, charging not possible
- Once charging is complete, disconnect the scooter from the charger and unplug the charger from the mains.

#### 6.3.2. Replacing the batteries

Be very careful when working on the batteries. Do not use metal tools, do not touch the battery terminals and ensure that you follow the battery manufacturer's instructions.


Caution! This work should only be done by a trained technician. There is danger of fatal injury!

Only use batteries that are approved for this scooter (see Technical data).

Make sure that the scooter is switched off before starting to replace the batteries. (See also Section 4.7 - Disassembling for storage).

- The batteries are located underneath the seat of the scooter and are protected by the housing cover.
- Start by removing the housing cover (Fig 41).
- - Disconnect the plug connections carefully and without using force.
- Undo the screw on the mounting bracket.
- Now remove the batteries (Fig. 43).
- Lift the new batteries into the space provided and connect the contacts. Make sure to observe the correct colour-coding (red and black, Fig. 42).
- Reattach the housing cover and the seat.

When the batteries have reached the end of their life, i.e. cannot be sufficiently charged any more, they must not be disposed of with the household waste. When new batteries are purchased, the mobility dealer takes the old batteries and disposes of them correctly.

## 7. Maintenance tips for retailers

We recommend having a complete inspection of the scooter carried out regularly by a qualified mobility dealer. If the scooter develops any faults or defects, it must be taken to the mobility dealer immediately for repair. The information, documentation and circuit diagrams required for repair and maintenance shall be supplied to the mobility dealer on request. A training programme is available for mobility dealers. If you have questions, you can contact the B+B customer service by telephoning +49 (0)700/6000 7070.

## 8. Storage / Shipping

If the scooter is to be stored or shipped, all attachable and unsecured parts must be removed and packed individually in appropriate cartons. The individually packed parts can then be packaged together in a larger carton. We recommend keeping the original packaging and storing it with these purposes in mind so that it is available if required. Your scooter then has optimum protection against environmental influences during storage or transport.

## 9. Disposal

If your scooter is no longer being used and is to be disposed of, please contact your mobility dealer.

If you wish to dispose of the product yourself, contact local waste management companies to enquire about the waste disposal regulations that apply to your place of residence.

## 10. Use by others

The Centuro S7 is appropriate for use by others. For use by others, the scooter must first be thoroughly cleaned and disinfected according to the accompanying care and hygiene instructions.

The product shall get checked by a qualified mobility dealer for general condition, wear and shall get repaired if necessary.

## 11. Maintenance

A maintenance plan can be provided to you by our customer service.

## 12. Service life

If used as intended, and the safety, care and maintenance instructions are observed, the expected service life is up to 5 years. Beyond this time the scooter can be used continuously, provided it is in a safe condition.

## 13. Warranty

Guarantee items refer to all defects of the scooter which can demonstrably be referred to material defect or manufacturing defect. If there are complaints you have to provide a completed certificate of guarantee. We provide the statutory warranty of 2 years for our products. A warranty period of 1 year applies for replacement devices and repairs. For batteries, a warranty period of 6 year applies.

Please also observe the B+B care, warranty, hygiene and maintenance instructions. Your qualified mobility dealer will be happy to provide these to you.

**Warranty certificate**

**Model: Centuro S7**

Model number:\*

\* (must be inscribed by the mobility dealer)

**Dealer:**

Date and seal

(Please note our general business terms and conditions)

\*Information labels are located at seatadjustment rod

Non-observance of the operating instructions as well as improper use, maintenance jobs and technical changes or adaptations without the permission of the Bischoff & Bischoff company, lead to guarantee loss and the loss of product liability in general.

### Tip!

In spite of careful research and editing, incorrect information may have crept into our operating manuals. Rights to technical changes reserved. Illustrations may include accessories. The current version of each instruction manual is ready for download at the B+B website:

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**B+B France S.A.R.L.**